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BASIC CONCEPTS

Three Working Modes

The edger makes **three basic operations** corresponding to a screen called "mode" that you can identify by a **symbol** located on right side of the screen.

Mode	Screen	
Tracing		<div data-bbox="909 560 1500 638" style="border: 1px solid black; padding: 2px;">1-CAPTURA FORMA POR TRAZADOR</div> <p>Allow you to enter the following parameters:</p> <ul style="list-style-type: none"> • job number • frame type • frame side to be traced • frame PD or Bridge • tracing type • job destination • job origin
Lens Layout/ Blocking		<div data-bbox="861 974 1508 1052" style="border: 1px solid black; padding: 2px;">2-CENTRADO Y BLOQUEO DE LENTES</div> <ul style="list-style-type: none"> • lens to be laid out • vision type • horizontal decentration • height decentration • segment width for near vision • lay-out/blocking mode • type of decentration
Edging		<div data-bbox="853 1299 1300 1377" style="border: 1px solid black; padding: 2px;">3-RECORTADO DE LENTES</div> <ul style="list-style-type: none"> • lens to be edged • lens type • lens finishing type • bevel type • lens final size • edging pressure • polishing

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Keys: Description

The following table describes each key of the

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Key	Function/description
	Enter numerical values.
	Confirm the entered numerical value.
	Decrement the selected numerical value.
	Increment the selected numerical value.
	Move the lens left during roughing and/or rimless.
	Move the lens left during roughing and/or rimless.
	Tighten/loosen the lens clamp. <i>Note:</i> This key exists also to the left of the edging station, for left-handed people.
	Maintenance function. Refer to Annex A, paragraph <i>Emptying the main tank</i> .
	Correct error while in entry process.

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
CREATING A JOB
Tracing Mode
Functions : illustration

Functions : description

The following table describes each function with regard to the tracing screen of the Silver.

Function	Selection	Description
Fct0		Tracing Mode
Fct1		Metallic frame
		Plastic frame
		Pattern or demonstration lens tracing
Fct2		Tracing the right side of a metallic frame
		Tracing the left side of a metallic frame
		Tracing the right side of a plastic frame
		Tracing the left side of a plastic frame
		Tracing the right side of a pattern
		Tracing the left side of a pattern
Fct3		Distance between centers (DBC)
Fct5		Automatic tracing
		Manual tracing: placing the stylus tip into the groove manually – slow tracing

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LAYING OUT AND BLOCKING
Lay-out/Blocking Mode
Functions : illustration

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The following table describes each function of the Silver lay-out/blocking

Function	Selection	Description
Fct0		Lay-out/Blocking Mode
Fct1		Lay-out/Blocking of the right side
		Lay-out/Blocking of the left side
Fct2		Distance vision
		Near Vision (for multifocals)
Fct3		Horizontal movement from the boxing center in millimeters
		Monocular Pupil Distance
		Binocular Pupil Distance
Fct4		Vertical movement from the boxing center in millimeters: A vertical decentration height from the boxing center
		BOX height: From the pupil center to the lowest vertical point of the lens, at the bottom of the frame.
		MIX height: From the pupil center to the first vertical point located on the edge and in the lower part of the lens.
Fct5		Segment width

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EDGING A LENS
Edging Mode
Functions : illustration

Functions : description
The following table describes each function with regard to the operation of the Silver

Function	Selection	Description
Fct0		Edging Mode
Fct1		Edging the right side
		Edging the left side
Fct2		Glass lens
		Organic lens
		Polycarbonate lens
		Trivex lens (depends on the type of edging wheel mounted)
Fct3		Rimless Finishing
		Front Side Bevel: The bevel is placed to the front of the lens.
		1/3-2/3 Bevel: The bevel is placed at 1/3 from the front face and at 2/3 from the rear face of the lens.
		1/2 - 1/2 Bevel: The bevel is placed in the center of the lens.
		Rear Side Bevel: The bevel is placed to the rear side of the lens.
		Controlled bevel: the bevel is positioned by the operator (see paragraph <i>Controlling the bevel of the lens before finishing</i> below)

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Function	Selection	Description
		Automatic bevel; based on a 1/3 – 2/3 bevel, the bevel is placed as close as possible to the front surface of the lens. The lens size is adapted to the curvature of the frame for the best possible fit. The circumference of the lens is the same as the circumference of the frame.
Fct5		Retouch: final lens size
Fct7		Minimum Edging Pressure: to edge a thin lens
		Medium Edging Pressure: to edge a medium lens
		Maximum Edging Pressure: to edge a thick lens
Fct8		Edging cycle for fragile lenses: to edge a very thin lens
		Water repellent lenses: special program to edge water repellent lenses.
Fct9		Lens Polishing: only for organic lenses
		No polishing of the lens

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